

Bottle for ophthalmic squeeze dispenser

Bottle is designed to be compatible with Aptar's Ophthalmic Squeeze Dispenser (OSD). The bottle and dispenser together provide a packaging solution for preservative-free formulations



OVERVIEW



- Designed for preservative-free formulation.



- Excellent alternative to blow/fill/seal unit dose vials. One 10 ml container holds the equivalent of twenty, 0.5 ml blow/fill/seal vials.



- 10 ml polyethylene (PE) bottle, sterilized by gamma radiation.



- Easy to dispense and use by the patient; intuitive handling and good squeezability.



- Bottles available in white and translucent options.

MARKETS



Eye Drop



Pharmaceutical



Medical Device

FEATURES

OSD construction



For additional information, please visit:
tekni-plex.com/healthcare