


Single-dose vials

Pre-molded single dose vials provide a higher level of precision, functionality and quality as compared to blow/fill/seal technology



OVERVIEW



- **Ideal for preservative-free liquids.**
Single-dose containers are suitable for liquids, gels or creams.



- Single-dose compact packaging is ideal to **reduce weight and volume during transportation.**



- **Easy to open, use and reclose.**



- **Dosing precision; no wasted product.**



- **Available in different shapes and capacities,** sized from 0.3 ml to 10 ml.



- **Injection molding creates a smooth orifice** with no burrs, flashes or sharp edges, enabling safe dispensing near sensitive areas, such as eyes or mucous membranes.



- **Ink-jet coding, text embossing and/or labelling is available.** Vials are also available in colors and can be decorated by pad printing.



- **Manufacturing in ISO 7 and ISO 8 cleanrooms.**



- **Sterilization available with gamma, beta or EtO.**

MARKETS



Eye Drop



Vaccine &
Oral Drug



Diagnostic



Inhalation



Allergy
Treatment



Veterinary



Cosmetic



Pharmaceutical



Medical
Device



For additional information, please visit:
tekni-plex.com/healthcare