

Pentafill® A

Automatic machines for filling and sealing single-dose plastic containers


Automatic machines for filling and sealing single-dose plastic containers at different production speeds.



OVERVIEW

The Pentafill A automatic machines are available in five models: the A5, A10, A15, A25 or A50 providing different cycles per minute.

- Automatic loading
- Compact machine; suitable for small rooms
- Easy-to-use machine with quick changeover capability and easy maintenance
- Maximum flexibility and versatility enabling you to fill different plastic containers and capacities
- Suitable for pharmaceutical products, diagnostics, veterinary, cosmetics, nutritional supplements and other products in liquid, gel or cream form
- Different filling systems available to fill different products viscosities
- Pentafill A5 - A10 - A15 supported formats: 0.33 - 0.6 - 1 - 2 - 3 - 5 - 10 ml




Pentafill® A

Automatic machines for filling and sealing single-dose plastic containers

Pentafill A5 – 10 – 15 Technical Data

	A5	A10	A15
Capacity	Up to 5 strips/min Up to 25 vials/min	Up to 10 strips/min Up to 50 vials/min	Up to 15 strips/min Up to 75 vials/min
Voltage	220V - 50/60 Hz	220V - 50/60 Hz	220V - 50/60 Hz
Power consumption	3 Kw	3 Kw	3 Kw
Compressed air consumption	450 L/min	450 L/min	450 L/min
Pressure required	> 6 bar	> 6 bar	> 6 bar
Chiller: water consumption	2-3 L/min	2-3 L/min	2-3 L/min
Chiller: MIN pressure	2 bar	2 bar	2 bar
Chiller: water temperature	14/16° C	14/16° C	14/16° C
Weight	520 kg	540 kg	560 kg
Dimensions (L x W x H) (mm)	1668 x 1685 x 1740 mm	1668 x 1685 x 1740 mm	1668 x 1685 x 1740 mm

Pentafill A5 – 10 – 15 Layout



For additional information, please visit:
tekni-plex.com/healthcare