

## Pentafill® A50

Automatic machine for filling and sealing single-dose plastic containers

Automatic machine for filling and sealing single-dose plastic containers in a sterile or non-sterile environment.



### OVERVIEW

The Pentafill A50 automatic machine provides 50 cycles per minute production speed in a compact footprint for filling and sealing your plastic containers in 5-vial strips.

- Automatic loading
- Compact machine; suitable for small rooms
- Easy-to-use machine with quick changeover capability and easy maintenance
- Maximum flexibility and versatility enabling you to fill different plastic containers and capacities
- Suitable for pharmaceutical products, diagnostics, veterinary, cosmetics, nutritional supplements and other products in liquid, gel or cream form
- Different filling systems available to fill different products viscosities
- Vial sizes include: 0.33, 0.6, 1, 2, 3, 5 and 10 ml




# Pentafill® A50

Automatic machine for filling and sealing single-dose plastic containers

## Pentafill A50 Technical Data

	A50
Capacity	Up to 50 strips/min. Up to 250 vials/min
Voltage	380V 3Ph - 50/60 Hz
Power consumption	15 Kw
Compressed air consumption	750 L/min
Pressure required	> 6 bar
Chiller: water consumption	2 - 3 L/min
Chiller: MIN pressure	2 bar
Chiller: water temperature	14/16° C
Weight	1320 kg
Dimensions (L x W x H) (mm)	3200 x 1360 x 1940 mm

## Pentafill A50 Layout



For additional information, please visit:  
[tekni-plex.com/healthcare](https://tekni-plex.com/healthcare)