

Pentafill® M

Semi-automatic machine for filling and sealing single-dose plastic containers



Semi-automatic filling and sealing machine ideal for small to medium batches.



OVERVIEW

The Pentafill M machine is available in two models: the M5 or the M10, providing different cycles per minute.

- Compact machine appropriate for small production spaces
- Suitable for pharmaceutical products, diagnostics, veterinary, cosmetics, nutritional supplements and other products in liquid, gel or cream form
- Flexible and versatile; different plastic containers and capacities can be filled and sealed
- Easy-to-use machine with quick changeover capability and easy maintenance
- Manual loading
- Different filling systems available to fill different products viscosities
- Pentafill supported formats: 0.33 - 0.6 - 1 - 2 - 3 - 5 - 10 ml




Pentafill® M

Semi-automatic machine for filling and sealing single-dose plastic containers

Pentafill M Technical Data

	M5	M10
Capacity	Up to 5 strips/min Up to 25 vials/min	Up to 10 strips/min Up to 50 vials/min
Voltage	220V - 50/60 Hz	220V - 50/60 Hz
Power consumption	3 Kw	3 Kw
Compressed air consumption	350 L/min	350 L/min
Pressure required	> 6 bar	> 6 bar
Chiller: water consumption	2- 3 L/min	2- 3 L/min
Chiller: MIN pressure	2 bar	2 bar
Chiller: water temperature	14/16° C	14/16° C
Weight	350 kg	370 kg
Dimensions (L x W x H) (mm)	1246 x 960 x 1740 mm	1246 x 960 x 1740 mm

Pentafill M Layout



For additional information, please visit:
tekni-plex.com/healthcare