


Pre-filled applicator

Pre-filled applicator for vaginal cream with tamper evident cap to ensure hygienic and safe dispensing



OVERVIEW



- **Capacity up to 6,5 ml.** Easy to use and ready for use. Dispenser is pre-filled, eliminating the need to fill the tube before use.



- **Accurate dispensing:** The single-use, pre-filled applicator enables accurate dosing.



- **The outer layer protects the cannula to ensure hygienic dispensing.**



- **Injection molding produces a smooth surface** that ensures no burrs, flashes or sharp edges for safe dispensing.



- **The tamper evident seal ensures inviolability of the content.**



- **Outer layer eliminates the need for an overwrap or a pouch.**

MARKETS



Pharmaceutical



Vaginal




Medical Device

FEATURES

How it works:

1. After removing the tamper-evident seal, peel with dispenser and roll the left part (3 petals) on his own axis to use it as a plunger.
2. Introduce the plunger and allow the product to be released from the container.



For additional information, please visit:
tekni-plex.com/healthcare