


Rectal applicator

Rectal applicator features a cannula with seven holes to ensure the correct release of the product in the designated area



OVERVIEW



- **Our rectal applicator conforms to European health, safety and environmental protection standards and can include CE marking, if required.** Our quality assurance department supports customers with the necessary documentation for the marking.



- **Safe dispensing:** The applicator enables accurate dosing.



- **Applicator is easily attached;** the M9 thread enables easy mounting of the applicator to the top of the tube.



- **The cover cap of the applicator ensures hygienic use.**

MARKETS



Pharmaceutical




Rectal



Medical Device

FEATURES

Cannula has seven holes to ensure the proper release of the product.



For additional information, please visit:
tekni-plex.com/healthcare