


Vaginal applicator

Injection molded vaginal applicators feature a smooth surface with no burrs, no flashes or sharp edges for user comfort



OVERVIEW



- **CE marked applicator.** If required, our vaginal applicator conforms with European health, safety and environmental protection standards and can include CE marking. Our quality assurance department supports customers with the necessary documentation for the marking.



- **Capacity from 1 ml to 5,7 ml.**



- **Dispenser box contains 1 to 16 pieces.**



- **Suitable for cream and gel formulation.**



- **Ergonomic design** for a practical and comfortable use.



- **Injection molding provides a smooth surface** with no burrs, flashes or sharp edges.

MARKETS



Pharmaceutical



Vaginal



Medical Device



For additional information, please visit:
tekni-plex.com/healthcare