

Silicone Extrusion Medical Tubing

High-quality, silicone extrusion medical tubing available globally from TekniPlex Healthcare's manufacturing facility in China.



OVERVIEW

Silicone is inert, biocompatible and an excellent choice for pump, and other applications due to its ability to consistently hold tight tolerances. Backed by world-class engineering and technical support, TekniPlex Healthcare's manufacturing facility in China is able to offer the global medical device marketplace high quality silicone extrusion medical tubing with short lead times. Available in platinum or peroxide silicone tubing along with other highly engineered silicone formulations.

Why consider TekniPlex Healthcare silicone extrusion tubing?

- New supply source
- Inert, biocompatible
- Tight tolerances
- Strictly regulated and monitored raw material supply
- Excellent mechanical and electrical properties
- Gas permeable
- Steam, radiation and EtO serializable
- Good chemical, fuel and oil resistance
- Global manufacturing footprint

Target Applications

- Catheters
- Feeding tubes
- Drug delivery
- Peristaltic pumps
- Medical implantable applications
- Bio pharmaceutical
- Automotive applications
- Industrial (electrical)



For additional information, please visit:
tekni-plex.com/healthcare